Mobi Mesh

Instruction Manual

Care for a Healthy Life

CONTENTS

English	1
Español	22
Português	44

MODEL NO:9R-001
PLEASE READ ALL INSTRUCTIONS BEFORE USE.

Dear Customer,

This manual contains important information you must know about the product. Please read it thoroughly before use. Please always keep this manual with you. The Nebulizer's flow rate will be different by different medicine characteristic; Besides, when medicine is in low temperature that flow rate will be decrease.

INTENDED USE

Medical purpose Intended user

This product is intended to be used for inhaling medication for respiratory disorders.

Legally certified medical experts, such as doctor, nurse and therapist, or healthcare personnel or patient under the guidance of qualified medical experts.

Environment

This product is intended for use in a medical facility, such as Hospital, clinic and doctor's office, in a room of general household, and in open-air environment with a roof.

Warranty

The Mobi mesh is two-year limited warranted (with exception of mouthpiece, medical cup, mask) to be free from defective workmanship from date of purchase. Any defective part(s)

will be repaired or replaced at APEX Medical's option if the unit has not been tampered with or used improperly during that period. Make certain that any malfunction is not due to inadequate cleaning or failure to follow the instructions. If repair is necessary, contact your authorized service personnel for instructions.

NOTE-Be sure to retain a dated proof of purchase document to verify unit is within 2-year warranty period.

NOTE-This warranty does not cover providing a loaner unit or compensating for costs incurred in rental while said unit is under repair.

Two years shelf life and six month expected service life is recommended of the single patient use accessories. Precautions Warnings described in the instruction manual should be observed.

- The mask system does not contain DEHP .

IMPORTANT SAFETY INSTRUCTIONS

READ THIS BEFORE USING

When using an electronic appliance, basic precautions should always be followed, including the following.

WARNING:

1. Close supervision is necessary when this appliance is used by, on or near children, handicapped person or invalids. Besides, the dosage and frequency of the use also has to be specified by the prescribing physician.
2. Do not operate under the blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
3. Use this appliance only for its intended use as described in this manual.
4. Do not use cotton swabs or other object to touch the mesh; otherwise, product will not operate.
5. Do not use different type batteries in the same time.
6. Do not carry or store the product which is still remaining medicine in medical cup.
7. For the type, dose, and regimen of the medication, be sure to follow the instructions of a doctor.

8. Do not immerse them in water.
9. Do not use disinfection solution containing sodium hypochlorite.
10. Do not inhale by using water in the medical cup.
11. Do not drop or apply strong shock to the Main unit.
12. Do not use any attachments or accessories not recommended by the manufacturer.
13. Discard the waste batteries and device according to local law.
14. Do not use the device on top of or close to other electronic device. FEATURE:
15. This product utilizes piezoelectric vibration to turn liquid drug into fine mist which is then sent to the respiratory tract and oral cavity by a venting system so that therapeutic result can be achieved.
16. Piezoelectric vibration which produces 100,000 ~120,000 vibrations per second.
17. Compact and easily to carry.
18. After 20 minutes continuously operation or when the medical cup is no medicine , the appliance will automatic shutdown. However, the no-medication detector has the possibility not function if that mesh blocked by medication.
19. Built-in low-battery indicator.
20. No matter use the product in different angle, that still can operate normally in short term. When put the product upside down, medicine cannot contact to mesh, the product still can operate normally in few seconds (it can be different if different medicine is in using).

NOTICE FOR SAFETY OPERATION

To ensure safety, please do follow all rules below:

1. Before using

- Use this appliance in the correct manners outlined in this manual.
- This appliance had to operate within the temperature of 10 C ~ 40 C, outside the temperature range was not suitable to operate the product.
 1. The medicine treatment should according to doctor's advice.
 2. If the product operates improperly, please discontinue to use it and send it for repairing immediately or consult your dealer.
 3. If the medical cup is not functioning very well (such as very low flow rate....etc.), before using, please replace it.
- When this appliance has not been used for a while, if any part (such as medical cup, mouthpiece or mask) seems dirty, please proceed with cleaning procedure described in the cleaning and maintenance procedure.
- Do not add medicine over the indicated line marked on medical cup (8c.c marking).
- Do not disassemble or make any modification on this appliance.

2. When using

- Children or disable patients should be assisted by adults or medical personnel when using this appliance.
- Do not use this appliance under a high humid environment (>85% RH)

- Do not use this product if the medical cup without any medicine.
- Keep the appliance in horizontal position and steady during operation to avoid spillage.
- Do not open the cover of medical cup during operation.
- Keep eyes away from the mist output.
- Treatment period please follow the doctor's advice.
- Different patients should use individual clean medical cups and mouth pieces to avoid cross infection.
- Some medicine do not nebulizes well, please consult with your doctor.
- This appliance has a low-battery indicator. When the indicator light is green, the power is good for working. When the power light turns to orange, the battery is become low, but you still can useproduct, however you may need to replace battery soon. If the power light is turns red it means voltage is no power to operate, please replace the battery immediately.
- During using this appliance, if any irregularity or discomfort is observed, please discontinue use this appliance.
- This product has a built-in timer which shut off automatically in 20 minutes. If it cannot be shut off automatically, discontinue using and consult your dealer for help.


- If product detects no medicine in the medical cup, the LED indicator will turn red and flash three times then shut off automatically.
- When using, you can use the product in any angles. But, make sure medicine has good contact with mesh. Please clean all parts after each usage and do follow instruction in the manual. Do not immerse this appliance in water when cleaning to avoid damaging any electrical parts.

WARNING: The turbo mode is for self clean purpose only, it can not be used in normal daily treatment, please consult your local distributor for further instruction.


3. After using

- Please clean all parts after each usage and do follow instruction in the manual.
- Do not immerse this appliance in water when cleaning to avoid damaging any electrical parts.
- Do not store this appliance in direct sunlight, outside our suggestion storage temperature and humidity environment.
- Keep this product from the reach of children.

COMPONENTES


1. Power button: Press to turn on or off the appliance.
2. Indicator Light: Indicate the power condition. (Green for normal/Orange for low/Red for no power). If there is no medicine inside, indicator light will turn red and flash three times before turn off automatically. However, the no-medication detector has the possibility of no function if the mesh is blocked by medication.
3. Medical Cup : Contain the medicine. (Suggestion:Replace a new medical cup per year)
4. Push button : Push here to release the medical cup from the main body.
5. Medical cup cover : Open here to pour the medicine.
6. Mouthpiece or Mask.
7. Mesh
8. Battery cover : Open here to replace batteries.

HOW TO USE

1. Make sure that all parts such as mouthpiece or mask, medical cup are clean, or clean the above three parts before using.
2. Open the battery cover and put in batteries. (The battery life varied by using different brands of battery)


1. Remove batteries if the unit is not in use for long periods. Failure to do so could result in damage due to battery leakage.
2. Discard the waste batteries according to local law.
3. Open the medical cup cover and pour the medicine into medical cup and close the cover completely.
4. Install mask or mouthpiece properly and start using the appliance.

- When the power button is pressed and the appliance starts working, the green light will on (sometimes, in orange and which depend on battery power remained) and the built-in timer will start to countdown.
- If the indicator turns orange, it shows the battery is low power. If the appliance is activated and the power indicator turns to red, it means the batteries are no power and the appliance cannot operate properly. Please replace with new battery immediately.
- If there is no medicine inside the medical cup, indicator light will turn red and flash three times before appliance turning off automatically.

CLEAN AND MAINTENANCE NOTICE

Before and after each using, please clean the following parts:

1. Remove the medical cup and mouthpiece or mask.

Discard all remaining medicine in the cup.

- Open the medical cap and discard all remaining medicine.
- Pour in some hot water and close the cover completely.
- Press the power button and operate 1~2minutes to discard all remaining medicine completely.

Ensure medicine be cleaned completely after each using, otherwise, mesh can be blocked.

2. Clean the surface of main body of appliance with cotton swabs or a piece of fine cotton fabrics.

Please clean the electrode in the main body and medical cup, ensure normally electrical conduction.

Please do not use easily volatility liquid to clean the electrode, ex.benzene, thinner.

3. Clean the medical cup, mouthpiece and mask with clean water and store at a safety place after cleaning the surface with ethanol.

Suggestion: Use a new medical cup for different illness or patients to avoid cross infection.

DISINFECTION

Medical cup can be disinfected by boiling in hot water (90°C):1 min, Validated number of cycles are 20. After 1 minute boiling, shake off excessive water from the cup, and let dry on a clean, lint free(paper) towel. Main unit, mask, mouth piece cannot be boiled.

TROUBLE SHOOTING

After checking the followings, if any problem still persists, please return the appliance to your dealer for servicing.

Q1. Press the power button, but no operation


1. Check if the battery voltage is sufficient? When the appliance's red indicator light is on, the battery is no power. Please replace the battery immediately.
2. Discard the waste batteries according to local law.
 - Check the battery is installed correctly?

Open battery cover and check if the battery electrode is rusted or not.

Q2. Cannot nebulize or not enough mist output

1. Check if the medical cup is filled with appropriate amount of medicine?
2. The mesh in medical cup blocked by remaining medicine. Cleaning medical cup, if still can't solve problem, please replace a new medical cup.
3. The mesh surface(mesh output)has covered by water or medication. Remove water or medication from mesh surface.
4. Medical cup is not installing properly. Release the medical cup from main body and push it back again.

Suggestion : Replace a new medical cup per year.

THE ACCESSORIES

Item	Mode no.	Size	Single patient reusable
Child Mask	9R-001MS	87.22 x 64.35 x 57 mm	YES
Adult Mask	9R-001ML	124 x 77 x 68.5 mm	YES
Mouth Piece	9R-001MP	43 x 27.5 mm	YES
Medical Cup (Gray)	9R-001MCG	54 x 43.99 x 38.50 mm	YES
Adaptor (USA),(not included)	9R-AU	90 x 76 x53 mm	N/A
Adaptor (EU),(not included)	9R-AE	90 x 76 x53 mm	N/A

Note: In order to provide the optimal performance and avoid device damage, please use our Genuine AC adaptor. * Our adaptor polarity is represented by the following symbol.

ELECTROMAGNETIC COMPATIBILITY

Due to numbers of electronic devices are used in your near environment and when medical device is in use which may be susceptible to electromagnetic interference from these electronic devices. Electromagnetic interference may cause incorrect operation of the medical device and could generate unsafe situation. In order to prevent unsafe situations, the device is designed and tested following EN60601-1-2 standard which defines the levels to immunity to electromagnetic interference as well as maximum levels of electromagnetic emissions for medical device. Even though the electromagnetic emissions of the device are extremely low, special precaution still need to be observed:

- It is recommended that you do not place the device on top of or close to other devices (such as mobile phone). Should you notice any interference with other electrical devices, keep the device away from those electronic devices.

SYMBOLS DESCRIPTION

Class II

Attention, should read the instructions.

IP22 Protected against solid foreign objects of 12,5 mm and greater; Protection against vertically falling water drops .when the enclosure is tilted at an angle up to 15° from its normal position.

BF symbol, which indicates that this product complies with degree of protection against electric shock for type BF equipment.

Direct current

Read the documentation symbol

Batch Symbol (yyymmddxxx)
yyy:ROC calender year,mm:month,d:day,xxx:lot no.

Serial number Symbol

Manufacturer

Authorised representative in the European community.

Refer to instruction manual

CE mark and notified body registration number, the Annex II, excluding Section 4 of EC directive 93/42/EEC have been met.

Reference number

Disposal of Electrical & Electronic Equipment (WEEE):

This product should be handed over to a suitable collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city council, household waste disposal service or the retail store where you purchased this product.


Temperature limitation/temperature range

Humidity limitation/Humidity range

TECHNICAL DATA

Model Number	9R-001
Power Supply	1.5VX 2, AA size alkaline battery AC adaptor (Input: 100~240VAC, 0.3-0.15A, 50-60Hz/ Output: 3VDC, 1.33A)
Power Consumption	< 2W
Viberating Frequency	103~123 kHz \pm 1K
Battery life	Up to 1.5 hr if use continously 2 weeks if daily use for 6 min treatment cycle
Battery Indication	Full (green), Low (orange), No power (red)
Flow Rate	0.35 \pm 0.1 mL/min
Particle Size	<5 μ m
AUTO-OFF Time	20 min
Nebulizer AUTO-OFF	No medication and no power
Medical Cup Volume	8mL (max.)
Weight	138g (exclude battery)
Dimension	44.5mm (L) x 53.4mm (W) x 121.8mm (H)
Medical Cup Volum	8mL (max.)
Accessory	Medical cup, mask (adalt & child), mouthpiece
Operation environment	10~40°C, 30~80%RH
Storage/Transportation	-20~70°C, \leq 85%RH
Sound	Noise level(at 1m distance) 50 dB

Performance may vary with drug such as suspensions or high viscosity .

More information can be obtained from the relevant medicine manufacturer.

Vapor Expert(PY-001) Portable Nebulizer - Instruction Manual

Recommended separation distances between portable and mobile RF communications equipment and the ME equipment

The ME equipment is intended for use in an electromagnetic in which radiated RF disturbances are controlled. The customer or the user of the ME equipment can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ME equipment as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d=1,2\sqrt{P}$	80 MHz to 800 MHz $d=1,2\sqrt{P}$	800 MHz to 2.5 GHz $d=2,3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

The ME declaration – electromagnetic immunity

The ME system is intended for use in the electromagnetic environment specified below.

The customer or the user of the ME system should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	Portable and mobile RF communications equipment should be used no closer to any part of the EQUIPMENT or SYSTEM including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Interference may occur in the vicinity of equipment marked with the following symbol. 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	

Declaration – electromagnetic immunity			
The ME system is intended for use in the electromagnetic environment specified below. The customer or the user of the ME system should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % UT (> 95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles < 5 % UT (> 95 % dip in UT) for 5 sec	< 5 % UT (> 95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles < 5 % UT (> 95 % dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EQUIPMENT or SYSTEM requires continued operation during power mains interruptions, it is recommended that the EQUIPMENT or SYSTEM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Declaration – electromagnetic emissions

The ME system is intended for use in the electromagnetic environment specified below.

The customer or the user of ME should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
Emissions CISPR 14-1	Passed	The ME system must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
Harmonic emissions IEC 61000-3-2	Passed	The ME system must be used only in a shielded location with a minimum RF shielding effectiveness and, for each cable that exits the shielded location, a minimum RF filter attenuation of 80 dB from 10 MHz to 20 MHz, 100 dB from 20 MHz to 80 MHz and 80 dB from 80 MHz to 100 MHz. (The minimum at 20 MHz is 100 dB and the minimum at 80 MHz is 80 dB.) The ME system, when installed in such a shielded location, is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Passed	

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Print-2015 All rights reserved
Apex Medical Corp.