**ADVANCE 1000**

**ADVANCE 2000**

*USER'S MANUAL*  
**SEQUENTIAL COMPRESSION THERAPY**

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*BEDIENERHANDBUCH*  
**SEQUENTIELLE PRESSOTHERAPIE**

**Model No.: AP-205000**

Please read the instruction manual before use.

# ***USER'S MANUAL***

English

User's manual .....P. 1

Español

Manual de Instrucciones .....P.8

Português

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Français

Instructions d'utilisation .....P.25

German

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## COMPRESSION THERAPY

The sequential compression therapy system has designed for lymphatic drainage. The sections in the sleeves are made of partially overlap so as to eliminate any gap between one another and reduce the possibility of stagnant fluid. The sections in the sleeves inflate and deflate sequentially. Compression starts in the distal chamber and moves on to the proximal one, so favoring fluid drainage. The pressure can be adjusted by means of a knob situated on the front of the pump.

### **Two systems are available:**

*ADVANCE 1000* with a 3-outlet pump for full leg or full arm sleeves.

*ADVANCE 2000* with a 2-outlet pump for half leg or half arm sleeves.

Both systems are lightweight and easy to operate.

## SYMBOL Definitions

Authorized representative in the European community.

Manufacturer

BF symbol, which indicated this product, is according to the degree of protecting against electric shock for type BF equipment.

**CAUTION** should be used to highlight the fact that there are specific warnings or precautions associated with the device.

Consult instructions for use

Attention – Observe proper Disposal of Electrical & Electronic Equipment (WEEE):

This product should be handed over to an appropriate collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

## ***EMC Warning Statement***

This equipment has been tested and found to comply with the limits for medical devices to the EN 60601-1-2:2007. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.

## ***General Information***

The system is a high quality and affordable mattress system suitable for treatment and prevention of pressure ulcers.

The system has been tested and successfully approved to the following standards:

EN 60601-1  
EN 60601-1-2

## **INSTALLATION**

### **BEFORE USE**

1. Before use make sure that the feeding cable is not damaged or crushed.
2. Supply power to the pump and turn it on without connecting the tubes. Check that air comes out of the outlets of the tubes (caution: the air comes out alternately from the outlets which can be in 2 or 3 according to the model).
3. Check that the connecting tubes are not damaged or crushed.
4. If necessary, replace the damaged parts only with original ones.

## **INSTRUCTIONS**

### **Turn the appliance off**

1. Place the pump on a flat surface or hang it where you can by means of the hooks.
2. Connect the tubes to the pump and to the arm/leg sleeve.
3. Check that the tube connecting points on the sleeve is not damaged.
4. Turn the pressure-adjusting knob to the minimum.
5. Put the sleeve on and fasten the zipper completely.
6. Switch on pump on. The green light indicates that the pump is working. Do not put on or take off the sleeve while the pump is working.
7. Adjust pressure turning the knob on the front of the appliance clockwise to the value indicated by your physician according to the following chart.

Knob	1	2	3	4	5	6	7	8
MmHg	30	38	47	55	64	72	81	90

The chart indicates the millimeters of mercury as units of measure and pressure and the corresponding numbers on the pressure-adjusting knob.

### **ATTENTION:**

**Do not open the zipper with pressure in the arm-leg sleeves. Before opening or closing the zipper make sure that the system (arm-leg sleeves) are deflated.**

## **TROUBLESHOOTING**

### **IN CASE OF DYSFUNCTIONING CHECK THE FOLLOWING POINTS.**

1. Check that the pump is connected to mains and that the green switch is on and lit.
2. Check that the tubes are properly connected to the sleeve and to the pump.
3. The sleeve does not inflate immediately. You have to wait for a few minutes before noticing the cyclical working of the pump.

**CONTRAINDICATIONS**

COMPRESSION THERAPY IS CONTRAINDICATED FOR PEOPLE SUFFERING FROM FOLLOWING PATHOLOGIES: Cardiopathy, suspected or ascertained thrombosis, arteriosclerosis or ischemia, dermatitis, gangrene, local infections.

**CLEANING INSTRUCTIONS**

Before cleaning the appliance, make sure it is not connected to the electric network. Do not clean the appliance while you are using it. To clean the pump, only use water and deterging soap. Avoid using solvents. You can wash the sleeve with liquid detergents in water not exceeding 50° C. Do not scrub strongly and never immerge the tube connectors into water. Do not dry-clean or iron.

Sterilization can be carried out only with gas at a temperature never exceeding 50° C.

	<b>ADVANCE 1000</b>	<b>ADVANCE 2000</b>
Power:	230V/50Hz	230V/50Hz
Pressure:	30 - 80 mmHg	30 - 80 mmHg
Air output:	5 liters	5 liters
Cycle time:	3 - 6 min.	3 - 6 min.
Feeding cable:	2,5 m	2,5 m
Outlets:	3	2
Dimensions:	24 x 13 x 10 cm	24 x 13 x 10 cm

Type:	<b>Full leg</b>	<b>Full arm</b>
Material:	Nylon + PU	Nylon + PU
Dimensions:	84.7 x 32 cm	78 x 31.5 cm

Type:	<b>Half leg</b>	<b>Half arm</b>
Material:	Nylon + PU	Nylon + PU
Dimensions:	62,5 x 35 cm	52,5 x 27,5 cm

## APPENDIX A: EMC INFORMATION

### Guidance and Manufacturer's Declaration- Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
Harmonic emissions IEC61000-3-2	Class A	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network.
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	

### Guidance and Manufacturer's Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Electrostatic Discharge (ESD) IEC61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	±2kV for power supply line ±1kV for input/out line	±2kV for power supply line ±1kV for input/out line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0,5 cycle 40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0,5 cycle 40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	Mains power quality should be that of atypical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of atypical location in a typical commercial or hospital environment.
NOTE: $U_T$ is the a.c. mains voltage prior to the application of the test level			

### Guidance and Manufacturer’s Declaration - Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Conducted RF IEC 61000-4-6	3Vrms150 kHz to 80 MHz outside ISM bands <sup>a</sup>	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance  <math>d = 1.2\sqrt{P}</math> 150kHz to 80MHz</p> <p><math>d = 1.2\sqrt{P}</math> 150kHz to 80MHz  <math>d = 2.3\sqrt{P}</math> 80 MHz to 2.5G MHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).<sup>b</sup></p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>c</sup>, should be less than the compliance level in each frequency range<sup>d</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.

b) The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.

c) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land

mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

d) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

### Recommended separation distances between portable and mobile RF communications equipment and this device:

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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Print-2018/All rights reserved