



Carry 185 Classic

patient hoist
designed for home care



sl-fastener
and adjustable handle



handset with magnet



spring lock



● MADE
● IN
● GERMANY

- gentle pick up and settle of the patient
- easy to use in narrow spaces and doors
- adjustable spreading for easy access to chairs and wheelchairs
- protects the back of the caregiver
- user-friendly due to low total weight
- lifting from floor level possible
- standard mast and chassis are removable without tools
- comfort-control with display
- soft-start / soft-stop, handling via control box possible

facts

appropriate slings:

- all slings suitable for 4 position points

article-no.:

- (euro pallets compliant packed) manual spreading: **8002**

optional:

- Carry 4-point-spreader-bar: **8060**
- Carry 4-point-spreader-bar basic: **8064**

standard:

- Carry 2-point-spreader-bar

material and service

material:

- powder coated steel
- ball-bearing swivel-casters

medical aid number: (Germany)
manual spreading: 22.40.01.0099

data

manual spreading

length	115 (45.2")
height	137 (54")
width	61 (24")
leg height	11 (4.3")
lift height	81-203 (31.8"-80")
rotating radius	132 (52")
spread	51-94 (20"-37")
weight	37 kg (81 lbs)

max. patient weight: **185 kg (410 lbs)**
further details on page: **61**

measure in cm