


Carry 185 Eco

patient hoist
designed for home care



sl-fastener
and adjustable handle



spring lock



Carry 2-point-spreader-bar



● MADE
● IN
● GERMANY

- gentle pick up and settle of the patient
- easy to use in narrow spaces and doors
- adjustable spreading for easy access to chairs and wheelchairs
- person recording from the ground possible
- foldable without tools
- underrun height front 110 mm

facts

appropriate slings:

- all slings suitable for 4 position points

article-no.:

- (euro pallets compliant packed)
manual spreading: **8000**

optional:

- Carry 4-point-spreader-bar: **8220**

standard:

- Carry 2-point-spreader-bar

material and service

material:

- powder coated steel
- slide-bearing swivel-castors

medical aid number: (Germany)
manual spreading: 22.40.01.0098

data

manual spreading

length	115 (45.2")
height	137 (54")
width	61 (24")
leg height	11 (4.3")
lift height	81-203 (31.8"-79.9")
rotating radius	132 (51.9")
spread	51-94 (20.8"-37")
weight	36 kg (80 lbs)

max. patient weight: **185 kg (410 lbs)**

measure in cm