



# Icon Floor

Minimising the risk of fall and potential harm

**APEX**  
Care for a Healthy Life

# Icon Floor

Minimising the risk of fall and potential harm

The electronically controlled Icon Floor bed is designed specifically with the occupant in mind, combining quality construction and functionality to assist with care, as well as reducing the risk to the occupant of falling and causing injury by using the low height functionality.

The stylish wooden bed has additional features including backrest, leg-lift and an overall bed height adjustable from 10cm to 70cm to further aid mobility and care.

Using a combination of these functions, the Icon Floor can also be placed into a comfortable profiling position for daily activities such as watching television or interacting with visitors making it the perfect profiling bed for multiple care environments including home, residential and nursing care.

Static foam and dynamic mattresses available.



## Key Features

- Minimum mattress base height: 10cm from the floor
- Easy to use handset
- 2 stage lowering to prevent entrapment
- Electric angle adjustment of the foot section
- Auto contour function
- Auto regression function
- Easy clean and maintenance
- Aesthetically designed head/foot boards for enhanced visibility
- Independently braked castors
- Quick and easy assembly
- 3 year warranty on metalwork and electrical actuator

## Specifications

Maximum user weight	: 180KG
Lowest platform height	: 10cm from floor
Maximum height position	: 70cm
Mattress platform	: 89cm- length 200cm
Extended	: 11KG
Castor diameter	: 3" (PVC) WH001
Backrest angle	: 73° Degrees
Calf section angle (when high section raised)	: 25° degrees
Auto regression function	
Trendelenburg/Reverse Trendelenburg	
Auto profile function	
Anti-entrapment function when lowering below 30cm	
Static foam and dynamic mattresses available	

### UK

Apex Medical Limited  
Unit 33, Great Western Business Park  
McKenzie Way, Worcester, WR4 9GN  
T: +44 (0)1905 774 695  
sales@apexmedicalcorp.co.uk  
www.apexmedicalcorp.co.uk



www.apexmedicalcorp.co.uk

Distributed by