



Pro Bario

Bringing Comfort to Bariatric Users

APEX
Care for a Healthy Life

Pro Bario

Bringing Comfort to Bariatric Users

The Pro-Bario is a very low bariatric care bed, with an adjustable height range of nominally 22cm to 62cm, minimising the risk of falls and enabling easier access to and from bed.

Two weight capacity options are available to suit user needs, and the bed's durable construction ensures maximum stability and support for the bariatric user. Pro-Bario also features an innovative anti-tilt system, preventing distortion and encouraging stability in the chassis over time and therefore prolonging lifespan.



Pro Bario
200cm (L) x
120cm (W)



Pro Bario+
206cm (L) x
127cm (W)



Pro Bario Max
200cm (L) x
120cm (W)

Key Features

- Easy assembly and Installation
- Electronically adjusted height from 22cm to 62cm
- Electronically adjusted Trendelenburg/Reverse function.
- Innovative anti-lateral tilt system for total stability
- 4 braked twin-wheeled castors
- Beech timber construction
- 4 section mattress support
- Powder-coated steel construction for durability and ease of cleaning

Specifications

Specification	Pro Bario	Pro Bario +	Pro Bario Max
Safe Working Load	260kg	260kg	320kg
Max User Weight	220kg	220kg	280kg
Platform Height	22-62cm	22-62cm	22-62cm
Mattress Base Width	120cm	127cm	120cm
Mattress Base Length	200cm	206cm	200cm
Backrest angle	0° to 70°	0° to 70°	0° to 70°
Leg	0° to 48°	0° to 48°	0° to 48°
Fowler Knee Break - 7 positions	4° to 26°	4° to 26°	4° to 26°
Trendelenburg / Reverse Trendelenburg			
Static Foam and Dynamic Mattresses available			
Actuators Ingress Protected IP54			

UK

Apex Medical Limited
Unit 33, Great Western Business Park
McKenzie Way, Worcester, WR4 9GN
T: +44 (0)1905 774 695
sales@apexmedicalcorp.co.uk
www.apexmedicalcorp.co.uk



www.apexmedicalcorp.co.uk

Distributed by