

# Relievo ProfiSoft Dual

## Pressure Redistributing Foam Mattress

- Suitable for patients in medium to high risk of pressure injuries
- A foam mattress replacement that can be used both sides for different patients as a cost-effective measure for institutions




Care for a Healthy Life

**PEX**

# Relievo ProfiSoft Dual

## Pressure Redistributing Foam Mattress



**Different Pressure Management on Two Sides**  
Two-density layer design in one mattress to meet with the diverse softness requirements for different patient weights



Maximum Patient Weight Up to 150 kg


Soft layer for patients between 0 to 70 kg

Firm layer for patients between 70 to 150 kg




**Greater Comfort and Microclimate Management**

The open cellular structure of Relievo ProfiSoft Dual offers a perfect ventilation beneath patient's body for better comfort level and microclimate management



**4-way Stretch PU**

Breathable and water impermeable PU cover for hygiene and easier maintenance



**High Density for Greater Comfort**  
50 kg/m<sup>3</sup> high density mid-layer foam for great comfort and high support to the patient.

### Specifications

Relievo ProfiSoft Dual		
Mattress	Type	15 cm-high foam mattress replacement
	Dimensions (cm)	200 x 90/100/120 x 15
	Top Cover Material	4-way Stretch PU
	Top Layer	45 kg/m <sup>3</sup>
	Mid Layer	50 kg/m <sup>3</sup>
	Base Layer	35 kg/m <sup>3</sup>
	Certificate	OEKO-TEX STANDARD 100
	Maximum Patient Weight	150 kg

*Specifications are subject to change for improvement without notice*

### Apex Medical Corp.

No. 9, Min Sheng St., Tu-Cheng,  
New Taipei City, 23679, Taiwan  
T +886 2 2268 3100  
F +886 2 2268 6525  
marketing@apexmedicalcorp.com

All rights reserved Specifications may change without notice.

www.apexmedicalcorp.com