



Sedens 400

Hybrid Seat Cushion

- Alternating cushion works with Apex dynamic mattress pumps to provide effective seating solution
- Lightweight and portable, perfectly to be used in chairs
- Hybrid design with alternating mode to redistribute pressure effectively
- Innovative air valve to control the amount of immersion for optimal pressure distribution



Care for a Healthy Life

.PEX

Sedens 400

Hybrid Seat Cushion



24 Hours Pressure Area Care
Alternating cushion works with Apex dynamic mattress pumps to provide effective seating solution



Hybrid Cells
6 air-cell-like one-piece pad with foam encased in to reduce risk of bottoming out



Positioning Cushion
Innovative air valve to control the amount of immersion for optimal pressure distribution

Non-Slip Base
Prevents from sliding out of position, giving peace of mind to the user

4-Way Stretch Cover
Fluid resistant, low shear and vapor permeable top cover protect skin from friction and moisture



Stretch Cover



Water Resistance



Flame Resistance



Biocompatibility



Vapour Permeable



Fungistatic



4-Way Stretch

Specifications

Cushion	Cushion Dimension	17.3" x 16.1" x 3.5" / 44 x 41 x 9 cm
	Cell Height	3.5"
	Weight	2.4 lb / 1.11 Kg
	Cover Material	4-Way Stretch
	Cell Material	Nylon / PU
	Base Material	Polyester PVC
	Arrangement	6 Hybrid Cells (Cellular bladders filled with foam)
	Flame Retardant Standards	EN597-1 & EN597-2
	Maximum Patient Weight	There is no weight limit, yet the cushion must be properly sized to the individual

Apex Medical Corp.

No. 9, Min Sheng St., Tu-Cheng,
New Taipei City, 23679, Taiwan
T +886 2 2268 3100
F +886 2 2268 6525
marketing@apexmedicalcorp.com