



Sedens 500

Seat Cushion with Battery Power

- Lightweight and portable, perfectly to be used in chairs or wheelchairs
- Solid foam base with medial thigh support for postural stability
- Lithium-ion battery of long duration that allows an autonomy up to 12 hours
- Visual and audible alarm enhances patients safety and security




Care for a Healthy Life

APEX

Sedens 500

Seat Cushion with Battery Power





Stretch Cover

Water Resistance


Flame Resistance

Biocompatibility

Vapour Permeable

Fungistatic

4-Way Stretch



A pump carry bag with adjustable strap can fit all sizes of wheelchair

Specifications

Pump	Dimension	9.8" x 4.7" x 2.2" / 25 x 12.2 x 5.7 cm		
	Weight	2.2 lb / 1.0 Kg		
	Case Material	Flame retardant ABS		
	Supply Voltage	12V DC, 1150 mA		
	Operating Cycle	10/15/20 min. selectable		
Cushion	Li-ion Battery Pack	3.78V Li-ion Battery x 2		
	Battery Operation Time	Approx. 12 hours		
	Battery Recharge Time	Approx. 4.5 hours		
Cushion	Dimension	16" x 16" x 4" / 41 x 42 x 10.2cm	17" x 17" x 4" / 43 x 43 x 10.2cm	18" x 18" x 4" / 46 x 46 x 10.2cm
	Weight	2.4 lb / 1.1 Kg		
	Material	Cover: 4-Way Stretch; Cell: Nylon/PU & PU Foam; Base: Nylon/PU		
	Maximum Patient Weight	260 lb / 120 kg		
	Flame Retardant Standards	EN597-1 & EN597-2		

Apex Medical Corp.

No. 9, Min Sheng St., Tu-Cheng,
 New Taipei City, 23679, Taiwan
 T +886 2 2268 3100
 F +886 2 2268 6525
marketing@apexmedicalcorp.com