



# Best Neb

## Comfortable Compressor Nebulizer

The Best Neb nebulizer is a high quality aerosol compressor that is designed to appeal to consumers, both in form and function. It is especially designed to avoid resembling medical equipment in order to make therapy more agreeable. This is especially important for pediatric care, where the majority of nebulizer users are children who suffer from the effects of poor air quality. Best Neb is the best choice in form, function and user friendliness.

# Best Neb

## Comfortable Compressor Nebulizer

### Features ▼

- Attractive appearance is less intimidating for children.
- Reusable inlet filter for simple, cost-effective maintenance.
- Anti vibration feet for stability on any stable work surface.
- Dispenses 3 cc's of medication in 10 minutes or less.



Convenient handle allows for easy portability and mobility.



Convenient storage compartment easily holds nebulizer accessories and tubing after every usage.



The clear lid keeps the contents visible for increased safety.



A hidden compartment allows the electrical cord to be stored away when not in use.

### Accessories Ordering Information:



**Mouthpiece set** : Anti-drool "T" mouthpiece with 6" reservoir tube and 7' supply tube.

**Mouthpiece set and 2 masks** : Anti-drool "T" mouthpiece with 6" reservoir tube, 7' supply tube and one each of adult & pediatric mask.

### Specification

Electrical Rating	120VAC, 60Hz, 1.2A	230VAC, 50Hz, 0.6A
(Note: Refer to the rating label on the product)	230VAC, 60Hz, 0.6A	
Max Pressure	≥35 psi	
Max Flow	≥14 lpm	
Operation Flow Rate	Approx. 5.5 lpm	Approx. 5.0 lpm
Nebulization Rate	0.2 ~ 0.3 ml / min	
Particle Size	0.5 ~ 5 microns	
Dimensions (W×D×H)	28.5 × 20.5 × 11.0 cm / 11.2" × 8.1" × 4.3"	
Weight	2.0 kg / 4.4 lb	

Distributed by

Apex Medical Corp.  
 No. 9, Min Sheng St., Tu-Cheng,  
 New Taipei City, 23679, Taiwan  
 T +886 2 2268 3100  
 F +886 2 2268 6525  
 marketing@apexmedicalcorp.com  
 www.apexmedicalcorp.com  
 Print-2015/04 / All rights reserved  
 Specifications may change without notice.