



SA210

Handheld Pulse Oximeter

If a handheld device at a very attractive price for monitoring patient's oxygen saturation in blood and heart rate with precision in all ranges from pediatric to adult, please look no further, APEX is the one.

SA210

Handheld Pulse Oximeter

Features ▼

- Touch-keypad for easy operation
- Backlit display for night or day
- High & low SpO2 and Pulse rate Program with Audible/visual Alarm
- Sound volume adjustment
- Pocket-size design. Compact and light for mobility
- 4 AA Alkaline batteries for power supply
- Automatic Power-off. Conserves battery life
- Data management software ; 3-mode versatility: Spot checking, memory, and PC-Link
- Orange pouch for the protection and mobility of the oximeter and probe



Specification

Scope of Measurement	SpO2 : 35% to 100%	Pulse rate: 30 - 250 bpm
Precision	SpO2 : 70% - 100% ± 2 digits(%)	35% - 69% : unspecified
Pulse rate	30 ~ 250±3 bpm	
Data Update period	less than 2 seconds	
Averaging	4 for SpO2 ; 8 for Pulse rate	
Weight	280g ; 0.33 ponds (excluding battery and probe)	
Size(HxWxD)	14.5 x 7.25 x 2.25 cm /5.70" x 2.85" x 0.89"	
Battery	4 "AA" alkaline /15 hrs continuous operation.	
Environmental Conditions	Operating Temperature: 32°F - 113°F (0°C - 45°C)	
Storage Temperature :	-4°F ~ 122°F (-20°C ~ 50°C)	
	Operating Humidity : 15% - 90%, noncondensing	
Standards	Type of Protection : Class II equipment / IEC60601-1-2,	
	Class B / EN60601-1, Type BF / ISO 9919: 2005	

Distributed by

Apex Medical Corp.
No. 9, Min Sheng St., Tu-Cheng,
New Taipei City, 23679, Taiwan
T +886 2 2268 3100
F +886 2 2268 6525
marketing@apexmedicalcorp.com
www.apexmedicalcorp.com
Print-2015/04 / All rights reserved
Specifications may change without notice.