



WiZARD 310/320

Nasal Mask / Full Face Mask
Sizing Gauge




WiZARD 310
Nasal Mask



L
M
S

WiZARD 320
Full Face Mask



1. Please print out this sizing gauge and ensure it is printed in scale.
2. Find your mask model and cut along dotted red line.
3. Place the fitting gauge on your nose refer to the left picture and determine appropriate cushion size.



This ruler is for scaling purposes.
Please print out this sizing gauge at 100%.

print as default (scale 1:1), do not fit to page